Vision
■ CLASSIC

TECHNICAL MANUAL

WOHLWEND AG

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Introduction

Since its introduction at the 1995 IDS in Cologne, Germany, Vision Classic has established itself as an easy to use ceramic with consistent shade accuracy. As a clinically proven ceramic, Vision Classic is now used in over 40 countries. Vision Classic is a special fine-structured ceramic developed for all alloys in a CTE range of 13,8–15,2 (25–600°C). Vision Classic also has a natural light transmission and reflexion. By using fluorescent and opalescent additional powders very individual and esthetic results are achieved.

Compatible with alloys having a CTE of 13.8 – 15.2 $\mu\text{m}/\text{mK}$ (25–600°C)

- Precious alloys
- Semi-precious
- Non-precious
- Non greening in combination with alloys containing silver

Wohlwend's unique manufacturing process gives an ideal grain size and pigment distribution. Benefits of this unique process include:

- Consistently accurate shade reproduction – V shades A1–D4
- Easy handling
- No slumping during build up
- Less shrinkage during firing

The vital shade reproduction and natural fluorescence of Vision porcelains correspond to natural teeth. This replication of nature can be obtained by either using a basic two layer buildup or by incorporating additional dentine and incisal modifiers into the restoration.

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Alloy Parameters

- Compatible with alloys
having a CTE of 13.8–15.2 (25–600°C).
- Ideally at 14.5 (25–600°C).
- If the CTE is higher 14.7 (25–600°C):
Use slow cooling cycle!
- If below 14.2 (25–600°C):
Remove restoration from lift as soon as furnace has fully opened.
- These parameters may vary between different alloys.
- In case of uncertainty we recommend to use Vision Bonding.

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Material Properties

We confirm that all porcelain powders stated on this list have been tested and conform to EN ISO 9693:2000 and 6872:1998.

Nome del prodotto Nombre del producto	Ceramica per metallo tipo 1 Porcelana para metal tipo 1				
Polveri Polvos		Resistenza alla flessione Resistencia a la flexión MPa	Solubilità Solubilidad µg/cm ³	Coefficiente di espansione termica Expansión térmica 25–500°C	Temperatura di transizione vetrosa Temperatura de transición vidriosa TG°C ± 10
Opaco in polvere Polvo opaco	Op	> 100 MPa	50 µg/cm ³	13,8	600°C
Opaco in pasta Opaco en pasta	PO	> 100 MPa	50 µg/cm ³	13,8	630°C
Dentina + Smalti Dentina + Incisal	D+S	> 80 MPa	25 µg/cm ³	13,3	590°C
Dentine opache Dentina opaca	OD	> 80 MPa	25 µg/cm ³	13,3	590°C
Intensivi per opachi Opaco intensivo	IO	> 100 MPa	50 µg/cm ³	13,8	600°C
Trasparenti Translúcido	T	> 75 MPa	25 µg/cm ³	13,3	590°C
Masse spalla Hombro	SP	> 80 MPa	25 µg/cm ³	13,8	590°C
Effetti incisali Efecto incisal	EF	> 80 MPa	25 µg/cm ³	13,3	590°C
Intensivi per dentine Dentina Intensiva	ID	> 80 MPa	25 µg/cm ³	13,3	590°C
Masse gengivali Gingival	ZF	> 80 MPa	25 µg/cm ³	13,3	590°C

Scala colori Color	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaco	O	O	O	O	O	O	O	O	O	O	O	O	O	O	O	O
Dentina Opaca	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD	OD
Dentina	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D	D
Smalto Incisal	I 58	I 58	I 59	I 59	I 60	I 57	I 59	I 59	I 59	I 60	I 59	I 59	I 60	I 60	I 59	I 59
Masse Transpa. Materiales translúcidos	TR: SCL (Superclear) – CL (Clear) – NT (Neutral) – OP (Opal) – R (Red) – BL (Blue) – GR (Grey)															

Modificatori di scala colori Modificadores de color	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Massa Spalla Material de hombro	SM neutral + BA				SM neutral + BB				SM neutral + BC				SM neutral + BD			
Massa eff. per Spalla Efecto para hombro	SM ivory, straw, yellow-gold, brown, yellow-grey, brown-grey															
Intensivi di croma per Opaco	O BA				O BB				O BC				O BD			
Opaco intensificador de chroma	10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensivi per caratterizzazioni per Op.	I O															
Opacos intensivos para caracterizaciones	W (White)				GB (Yellow)		OR (Orange)		OK (Ocher)		BR (Brown)		V (Violet)			
Intensivi di croma per Dentina	D BA				D BB				D BC				D BD			
Dentina intensificador de chroma	10%	20%	30%	35%	40%	10%	20%	30%	40%	10%	20%	30%	40%	20%	30%	40%
Intensivi per caratterizzazioni dentina	I D															
Dentina intensiva para caracterizaciones	EB (Ivory)		ST (Straw)		HY (Honeyyellow)			BR (Brown)		MG (Grey)		BL (Blue)		NEON		
Intensivi per Smalto	II															
Incisal intensivo para individualizaciones	7 8 9 10															
Masse effetto per mammelloni	EF ivory, straw, honeyyellow, khaki															
Efecto para mamelones																

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
Metal processing prior to porcelain application

1. Process alloys according to manufacturers instructions. Avoid sharp edges.

Unless specified otherwise, the alloy should be sandblasted with min. 110 μm aluminum oxide at 29 psi.

A uniform oxide layer is desired. In case of NP, the oxides are to be removed using min. 110 μm aluminum oxide. (Refer to manufacturers manual!). A clean and grease-free surface must be obtained.

2. Apply Vision bonding paste in a uniform and thin layer. Firing cycle according to the bonder manual.
3. The result after firing is also replacing a washfiring.

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First step: 1st Opaque

1. Powder / Liquid:
Mix the opaque to a creamy texture. Apply using brush or glass instrument with a 70% coverage. Firing temperature 1778°F.
2. After first firing.

Vis Opaco - Paste

Apply using brush or glass instrument. Firing temperature 1778°F. Can be reduced to degas temperature of alloys. Prolong hold time at final temp. 2-3 min. Secure well after usage to avoid a dry out of the material.

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Second step: 2nd Opaque

1. Apply the opaque with a full covering layer.
2. Egg shell appearance after firing.

Options


To achieve internal effects, characterizations can be made by using opaque-modifiers "IO".

As an example:

Layer IO-orange in the central Fossa or IO-brown in the cervical area. To reach a higher chroma from the depth, apply some Chroma opaque to the equatorial area or admix as per shadechart.

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Third step: Shoulder

Our shoulder powders in basic shades A, B, C, D and Neutral are giving you the chance to create a shoulder for every shade situation. Please follow the mixing ratio in the chart below.

Optionally you can reach a larger color palette by using effect shoulder masses.

Shoulder combination Basic Vision Classic:

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





Fourth step: Build up

1. In case of less space apply a thin layer of opacious dentin (occlusal, cervical, etc.).
2. Build up with dentin (shade that is given) in full contour.
3. Do a small cut back from incisal and approximal. Apply the borders and mamelons with a thin layer of transparent neutral.
4. Complete the anatomic form of the tooth with enamel according to the shade. A light oversize of the build up, brings a perfect result, because of the small shrinkage during the firing cycle.

Use some opacious dentine to fill up interdental space.

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Fifth step: Correction

1. After the first dentine firing.
2. If needed correct the form with Transpa Super Clear.
3. Proceed with second dentine firing.

Sixth step: Finishing

4. When the contact points and the occlusion are worked out, you have got two possibilities to glaze.

Glaze firing without glaze powder and subsequent a mechanical polish (special rubber wheels, brush, Dia glaze and so on). This is only recommended if the situation is well known.

Glaze firing with glaze powder. Mix it with glaze liquide and apply a thin layer. For individual characterization use Vision Classic stains.

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Multi layering technique

The use of Vision's extended range of modifiers in conjunction with simple layering techniques allow a natural looking restorations. These can be used straight or mixed into other materials.

Shade example: Orange neck, Body between A3 and A3.5, Incisal A2.



To raise the chroma mix 30% ID BA with 70% dentine A3.

Neck area Intensive body ID OR, Upper third opacious body OD A3.

Build up with dentin A3. Reduce dentin slightly more than usual.

Incisal S58 and transpa TR Opal are alternatly layered to achieve the shade but also getting natural opaleszenz.

TR Opal can be mixed with all transpas, depending on which level of translucency is to be achieved in correction.

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Multi layering technique

Mix the highly fluorescent ID Neon with the intensive incal IS 7, to achieve higher value and apply in a thin layer.

For mammelons our effect masses will be used very carefully and in a thin layer. Effectmasses (EF): Ivory, Straw, Honeyyellow, Khaki

Neon can also be mixed with appropriate stain and be applied onto wet surface. (What you see is what you get!)

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Classic Firing chart

Recommended firing parameters!
Temperatures can vary due to different furnace performances.

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* For the Vis-Opak pasteopaque extend the dry time for 4 minutes.

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16 Dentine powders A1–D4 20 g
2 Modellingliquid 25 ml

ASSORTMENT S3

Incisal: 57, 58, 59, 60 20 g; Intensiveincisal: 7, 8, 9, 10 20 g
Transpa: Clear, Superclear, Neutral, Opal, Red, Blue, Grey 20 g
Glaze 10 g; 1 Modellingliquid 25 ml; 1 Glazeliquld 25 ml

ASSORTMENT S4

16 Opaque-Dentine powders A1–D4 20 g
2 Modellingliquid 25 ml

ASSORTMENT S5

Intensivopaque: White, Yellow, Orange, Ocher, Brown, Violet 20 g
Intensive dentin: Ivory, Straw, Honeyyellow, Brown, Grey,
Blue, Neon 20 g; Gingiva: Softpink, Darkpink 20 g
1 Opaqueliquld 25 ml; 1 Modellingliquid 25 ml

ASSORTMENT S6

Shouldermasses: Ivory, Straw, Yellowgold, Brown, Yellowgrey,
Browngrey 20 g; Effectmasses: Ivory, Straw, Honeyyellow, Khaki 20 g
1 Shoulderliquld 25 ml; 1 Modellingliquid 25 ml

KIT S7

Opaco Base: O-BA , O-BB, O-BC, O-BD 20 g
 Spalla Base: SM-BA, SM-BB, SM-BC, SM-BD, SM-Neutral 20 g
 Dentina Base: D-BA, D-BB, D-BC, D-BD 20 g
 1 Liquido Opaco 25 ml; 1 Liquido Spalla 25 ml

KIT S8

Bleach Opaco: O-A0/BO 20 g; Bleach Dentina Opaca: OD-A0, OD-B0 20 g
 Bleach Dentina: D-A0, D-B0 20 g; Bleach Smalto: A0/BO 20 g,
 1 Liquido Modellazione 25 ml; 1 Liquido Opaco 25 ml

TEST KIT

O-A3 20g, OD-A3 20 g, D-A3 20 g, I-59 20 g, TR-Clear 20 g;
 1 Liquido Modellazione 25 ml; 1 Liquido Opaco 25 ml

STARTER KIT

O-A2 20 g, O-A3.5 20 g, O-B3 20 g, O-C3 20 g
 OD-A2 20 g, OD-A3.5 20 g, OD-B3 20 g, OD-C3 20 g
 D-A2 20 g, D-A3.5 20 g, D-B3 20 g, D-C3 20 g
 I-58 20 g, I-59 20 g, TR-Clear 20 g, TR-Opal 20 g
 1 Liquido Opaco 25 ml; 1 Liquido Modellazione 25 ml

STAIN KIT

01, 02, 03, 04, 05, 06, 07, 08, 09, 10, 11, 12, 13, 14, 15, Glasura
 Additional grey e violet; 1 Pennello; Liquido Glasura e Stain 25 ml

VIS-OPAK-KIT

16 Paste Opachi in Pasta A1–D4 4 g
 White, Ocher, Lila, Pink 3 g; 1 Pennello

ASSORTMENT S7

Chromaopaque: O-BA , O-BB, O-BC, O-BD 20 g
 Shoulderbase: SM-BA, SM-BB, SM-BC, SM-BD, SM-Neutral 20 g
 Chromadentin: D-BA, D-BB, D-BC, D-BD 20 g
 1 Opaqueliqid 25 ml; 1 Shoulderliqid 25 ml

ASSORTMENT S8

Bleachopaque: O-A0/BO 20 g; Bleachopaquedentin: OD-A0, OD-B0 20 g
 Bleachdentin: D-A0, D-B0 20 g; Bleachincisal: A0/BO 20 g,
 1 Modellingliqid 25 ml; 1 Opaqueliqid 25 ml

TESTKIT

O-A3 20g, OD-A3 20 g, D-A3 20 g, I-59 20 g, TR-Clear 20 g;
 1 Modellingliqid 25 ml; 1 Opaqueliqid 25 ml

STARTERKIT

O-A2 20 g, O-A3.5 20 g, O-B3 20 g, O-C3 20 g
 OD-A2 20 g, OD-A3.5 20 g, OD-B3 20 g, OD-C3 20 g
 D-A2 20 g, D-A3.5 20 g, D-B3 20 g, D-C3 20 g
 I-58 20 g, I-59 20 g, TR-Clear 20 g, TR-Opal 20 g
 1 Opaqueliqid 25 ml; 1 Modellingliqid 25 ml

STAINASSORTMENT

01, 02, 03, 04, 05, 06, 07, 08, 09, 10, 11, 12, 13, 14, 15, Glaze
 Additional grey and violet; 1 Brush; Glaze & Stainliqid 25 ml

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Accessory

- Porcelain crock
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Bonding

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**Wohlwend AG