
Vision  
■ ZIRKON

TECHNICAL MANUAL

WOHLWEND AG

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## Introduction

Innovations afford new materials! With its high strength and delicate frameworks, zirconium oxide is enjoying increased popularity. Wohlwend's Vision Zirkon provides a new standard for shade accuracy. Its exact CTE settings and excellent adhesion result in optimum performance on different zirconium oxide frameworks. This protects restorations from cracks, checks and fractures.

With our new press to zircon system we once more unite high aesthetics with constant quality and easily work engineering. The new development is compact and applicable for every dental technician. Valuable time is saved through the singular in fully contour principle and less needed for the rework. Hereby we achieve a higher efficiency with a high-quality result.

- Shade accuracy
- Superb natural esthetics
- Easy handling
- Great stackability / carving abilities
- Outstanding adhesion and excellent bonding power
- Exact CTE settings CTE 9.5

Press over Zircon

- pressable with or without ceramic shoulder
- natural color depth and translucency
- Stain color system in 4 V-Shades
- only one press liner
- Cut Back technique with Vision Zirkon possible

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## Material Properties

We confirm that all porcelain powders stated on this list have been tested and conform to EN ISO 9693:2000 and 6872:1998

Classe Clase	Solubilità Solubilidad	Resistenza alla flessione Resistencia a la flexión	Coefficiente di espansione termica CET (25/500°C)		TG  ± 10	Distribuzione della dimensione delle particelle D90/10 Distribución del tamaño de la partícula	Cristallinità rilevata con XRD Cristalinidad detectada por XRD
			max.   máx. 100 µg/cm <sup>2</sup> (± 5)	min.   mín. 50 MPa			
Liner	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita
Dentina, Dentina Opaca	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita
Smalti, Smalti intensivi Incisal, Incisal Intensivo	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita
Masse spalla Hombro	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	580 °C	45 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita
Trasparenti Translúcido	16 µg/cm <sup>2</sup>	> 70 MPa	9,5 ppm/K	9,5 ppm/K	550 °C	60 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita
Pigmento	–	–	9,5 ppm/K	9,5 ppm/K	550 °C	20 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita
Glasura Glaze	–	–	9,4 ppm/K	9,4 ppm/K	550 °C	20 µm	Vetro amorfo / Leucite Vidrio amorfo / leucita

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
To reach perfect results it is required to prepare the framework exactly and according to the manufacturers instructions.

### First step: Liner

A thin layer liner according to the shade will be put on the framework and fired, before starting second step.

### Second step: Shoulder

Three different shoulder powders are giving you the chance ...

- to create a shoulder for every shade situation.
- to lower the values mix in intensive Incisal white.
- to lower the chroma mix in transpa clear.

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



### Third step: Build up

1. In case of less space apply a thin layer of opacous dentin (occlusal, cervical, etc.).
2. Build up with dentin (shade that is given) in full contour.
3. Make a small cut back from incisal and approximal. Apply the borders and mamelons with a thin layer of transparent neutral.
4. Complete the anatomic form of the tooth with enamel according to the shade. A light oversize of the build up, brings a perfect result, because of the small shrinkage during the firing cycle.

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### Fourth step: Correction

1. After first dentin firing.
2. Work out the crown and correct the form with Transpa Clear.

To achieve opalescence in the incisal area, TR Opal can be mixed into all transpa materials.

Even on a simple build up, natural light transmission and reflection can be achieved.

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

## Fifth step: Finishing

1. When the contact points and the occlusion are worked out, you have got two possibilities to glaze.
2. Glaze firing without glaze powder and a subsequent mechanical polish (special rubber wheels, brush, Dia glaze and so on). This is only recommended if the situation is well known.


Glaze firing with glaze powder. Mix it with glaze liquid and apply a thin layer.

For individual characterization use Vision Zircon stains.

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## Press over Zircon

It is possible to press over all commercially Zirconia systems (green, sintered or HIPped). The grinded or sintered Zirconia framework has to be prepared with special grinders or rubberwheels to manufacturer's specifications. Here, the manufacturer's minimum strengths of copings and thickness of connections have to be followed.

### First step: Liner

1. Apply the special universal liner in a thin and watery consistence to the framework and fire as shown in the firing chart. This serves to reduce the brightness value of the zircon framework material, and to get a better bonding.

**Attention:** Different firing temperature.

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## Second step: Wax up

Use 100% ash free waxes only! Put the framework without pressure on the isolated model. It is possible to wax up in full contour, or to reach a higher esthetic work in a reduced way (cut back).

In such cases apply from our Vision Zircon incisals and transpa powders. There is no loss of dimension during pressing.

## Third step: Spruing

Use 6 or 8 gauge ash free sprues with a length of 4 to max. 7 mm. Sprues are connected to the coping in a straight line so as not to restrict the flow of the molten ingot. Sprues should be angled slightly outward with a 3 mm distance between copings.

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Wax weight 0.6–1.3 g = 2 Ingots

#### Optional:

Paper ring former or silicon ring former.

Remove strip from paper, form a ring to the appropriate mark. Slip onto base and adapt stabilizing ring. Pour special investment to just below the stabilizing ring.

Carefully remove the stabilizing ring and fit the ring former, twisting it carefully. Investment must leak from the opening. After bench set, remove paper / Silicon ring, base and ring former, then clean the bottom surface with a plaster knife. This must be absolutely flat, so that the ring stands perfectly upright in the pressing furnace.

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Pressare a 4,5–5 bar / 65–72 psi.

## Fifth step: Pressing


For full contour press technique use Transpa Dentin pellets. For cut back technique use Dentin pellets.

Observe pressparameters and combinationcharts.

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### Sixth step: Devesting

Mark the length using spare plunger. Cut along the line, using separating disk. Attention, use mask and dust suction. Break segments apart.

Carefully devest the objects by sandblasting. Use glass beads only.

Start with 60 psi pressure, as soon as objects appear turn pressure down to 35 psi proceed blasting carefully. **DO NOT USE ALUMINIUM OXIDE!**

Carefully separate the sprue from the objects with a sharp diamond disk.

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### Seventh step: Layering (cut back)


After working out the press channels, layer up with Vision Zircon. Complete the incisal part and fire according to the chart. With dentin you can layer supplemented portions. If necessary make a correction firing. Before glazing you can customize using Vision Zircon stains. Fire on ceramic pins. It is important to ensure that the ceramic shoulder has no contact with the pin. It can be glazed or polished.

### Eighth step: Staining technique

Individual nuances can now be added with Vision Zircon Stains.

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Temperatures can vary due to different furnace performances

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### KIT STAINS

White, Honey, Olive, Redbrown, Chestnut, Grey, Blue 3 g,  
Glaze, A, B, C, D, 3 g, 1 Liquido Glasura e Stain 25 ml

### STARTKIT

L-2 20 g, OD-White 20 g, D–A3 20 g, I-59 20 g, TR-Clear 20 g,  
1 Liquido Modellazione 25 ml

## Assortments

### ASSORTMENT S1

16 Dentine powder A1–D4 20 g  
2 Modellingliquid 25 ml

### ASSORTMENT S2

Incisal: 57, 58, 59, 60, 20 g  
Opalincisal: Neutral, Blue 20 g; Intesiveincisal: White, Yellow 20 g  
Transpa Clear 20 g, Glaze 20 g, Liner: 1 to 3 and Pressliner 20 g,  
1 Modellingliquid 25 ml, 1 Glaze & Stainliquid 25 ml

### ASSORTMENT S3

Transpa: Neutral, Yellow, Blue, Grey 20 g  
Opaquedentine: White, Honey 20 g,  
Gingival: Pink and Softpink 20 g  
Shouldermasses: Creme, Yellow, Redbrown 20 g  
1 Modellingliquid 25ml, 1 Shoulderliquid 25 ml

### ASSORTMENT S4

Pressables each 5 pc: A1, A2, A3, B1, B2, B3, C1, C2, D2, Bleach,  
TA1, TA2, TA3, TB1, TB2, TB3, TC1, TC2, TD2

### STAIN-ASSORTMENT

White, Honey, Olive, Redbrown, Chestnut, Grey, Blue 3 g,  
Glaze, A, B, C, D, 3 g, Glaze & Stainliquid 25 ml

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## Accessory

- Porcelain crock
- Brushes
- Mixing trays
- Ceramic isolation pen
- Firing paste
- Plunger investment
- Plunger Aluminium oxide
- Silicon ring former
- Investment

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**Wohlwend AG